




Amplicon-EZ Ordering Guide

1. Log in to your account and select the Amplicon-EZ order form.



2. Enter general information about your order. Required fields are marked with asterisks (*).

3. Click through the three steps of the table and enter your sample information.

STEP 1: BASIC INFORMATION

	Sample Name	Partial Adapter Sequences?	Sample Buffer	Size of Amplicon (bp)
	EXAMPLE	Yes	Tris/EB	325
1	Test_1	No	Tris/EB	400
2	Test_2	No	Tris/EB	400
3	Test_3	No	Tris/EB	400
4	Test_4	No	Tris/EB	400

STEP 2: CONCENTRATION

	Sample Name	Method of DNA Measurement	Normalized to 20 ng/μL?	DNA Concentration (ng/μL)	Sample Volume (μL)	Total DNA Amount (ng)	Gel Image?
	EXAMPLE	Qubit	Yes	20	25	500	Yes
1	Test_1	Qubit	Yes	20	25	500	Yes
2	Test_2	Qubit	Yes	20	25	500	Yes
3	Test_3	Qubit	Yes	20	25	500	Yes
4	Test_4	Qubit	Yes	20	25	500	Yes

STEP 3: DATA ANALYSIS

	Sample Name	Number of Targets	Reference Sequence	Forward Primer Sequence	Reverse Primer Sequence	Need Analysis?	Available Analysis
	EXAMPLE	One	ATGGTGAAGCAAG...	GTGTAAACGACGGCCA	TCACACAGGAAACAGCTAT	Yes	SNP/INDEL detection
1	Test_1	One	ATGAGTATCAACATTCCGGTGC GCCCTT			Yes	SNP/INDEL detection
2	Test_2	One	ATGAGTATCAACATTCCGGTGC GCCCTT			Yes	SNP/INDEL detection
3	Test_3	One	ATGAGTATCAACATTCCGGTGC GCCCTT			Yes	SNP/INDEL detection
4	Test_4	One	ATGAGTATCAACATTCCGGTGC GCCCTT			Yes	SNP/INDEL detection

4. Click **Save and Review**

5. Review the order details and click **Add to Cart**

6. Proceed to checkout. **Make sure to print out the order receipt and include it with your samples.**

